

# Action Gel Professional Cushion



The Professional provides exceptional pressure relief for users at risk of skin breakdown or with a history of pressure ulcers. A full 4.1cm thickness of Akton polymer provides effective weight distribution on a contoured or flat surface. The

Professional provides an additional 50% cushioning against pressure, shear and vibration. This is significant for those active clients using a power wheelchair who have increased concerns about bottoming out and the potential effects of increased temperatures (due to foam insulation properties). Custom sizes available upon request.

## PRODUCT BENEFITS:

- 41mm thick
- Suitable for very high risk clients, including medium risk as a preventative measure
- Not restricted by any weight limit
- Easy to care for
- Basic cover included

## AVAILABLE SIZES:

16" x 16"	3540
18" x 16"	3541
18" x 18"	3542
22" x 18"	35422
20" x 18"	3543
20" x 20"	3544

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au