

Deluxe Ultralight Alloy Wheelchair



The ActiveCare® Deluxe Ultralight Alloy Wheelchairs have been designed to accurately meet the needs of our clients. They are an ultralight weight wheelchair with the benefit of many additional features.

The improvement compared to previous designs (ranging from 25 - 30 kg gross weight) is considerable. Much input from clients and nurses has assisted in the development of these superior quality wheelchairs. The Deluxe Alloy wheelchair has all the features you need.

Features and Specifications



Flip back Armrests



Removable Armrests

BENEFITS & FEATURES:

- Flip back armrests
- Swingaway, detachable legrests
- Fold down backrest
- Quick release rear wheels
- 24" rear wheels with urethane tyres
- Aluminium frame
- Padded nylon upholstery
- Complies with AS3695:2008



Height Adjustable Armrests



Urethane Tyres



No Upholstery Screws

PRODUCT SPECIFICATIONS:

TECHNICAL DATA		
Product Code	4029	4629
Seat Width	400mm	460mm
Seat Height (from floor)	480mm	480mm
Seat Depth	400mm	400mm
Backrest Height	430mm	430mm
Armrest Height	230mm	230mm
Overall W x D x H	600 x 1060 x 960mm	650 x 1060 x 960mm
Fold Width	320mm	320mm
Wheel Diameter	24"	24"
Net / Gross Weight	11 / 14kg	11 / 14kg
SWL	125kg	125kg

